
Detecting Cancer Metastases on Gigapixel Pathology Images

Final Project of
Applied Deep Learning


Rui Bai (rb3454) Yichi Liu (yl4327)

Project Introduction

- Motivation
- Flow of the Project

Motivation

- Metastasis detection :
 - Detect whether the breast cancer has spread to nearby cells
 - Early diagnosis will help doctors to give treatment
- However, manually labelling the cell will be time-consuming.
- We designed an automatically Metastasis detection model on the Pathology image with CNN models



Cited: Detecting Cancer Metastases on Gigapixel Pathology Images, 2017

Flow of the Project


- Project Flow Introduction: introduce the overall steps as below
- Data Processing
 - Training & Validation & Testing
 - Patch extraction for training & validation data
- Model Architectures
 - Different transferred models
 - Different Scales
- Heatmap Construction
- Model Comparison
- Comparison of Results
- Final Prediction for 3 testing data

Data Processing


- Training & Validation & Testing
- Patch Extraction for Training & Validation data

Training & Validation & Testing

- Some slides have little cancer cell.
We don't take them into consideration
- Training and Validation Image:
Patches from 8 slides:
slide 016, 031, 064, 075, 078, 084, 094, 101
 - Training: 80% of patches
 - Validation: 20% of patches
- Testing Image:
Patches from 3 slides:
Slide 091, 096, 110
- Observation: **Imbalanced** dataset



Slide 012:
Little cancer cell.
Drop it.




Slide 035:
Almost no cancer cell.
Drop it.

Positive Patch Extraction for Training & Validation data

Randomly Get 200 Positive
299*229 patches:

- If the center point of the patches is not cancer: **Drop it**
- If the center point of the patches is cancer: **Save** the patch and label as **Positive**




Condition a:
No cancer at the center point.
Solution:
Drop it

Condition b:
The center point is cancer.
Solution:
Label is **POS**

Negative Patch Extraction for Training & Validation data

Randomly Get 200 Negative
299*229 patches:

- If the center point is not tissue (intensity > 0.8): **Drop it**
- If the center point is tissue:
 - If the center region (128*128) contains cancer: **Drop it**
 - If the center region doesn't contain cancer: **Save** the patch and label as **Negative**



Condition a:

No tissues at the center point.

Solution:

Drop it

Condition b:

The center region (128*128) contains cancer.

Solution:

Drop it

Condition c:

The center region (128*128) has no cancer.

Solution:

Label is **NEG**

Avoiding Memory Issues


- Don't Read the whole slide image
- Extract the mask and the slide image only with size 229*229
- Save the patches into folders




We could run our script with even level 0 without crashing in Colab
(without update to Colab Pro)

Negative Extraction for Training & Validation data

- Randomly Get 200 Positive patches:
 - Conditions:
 - If
- Negative patches:
 - Get tissue pixels
 - Randomly select a tissue pixels to capture a patch
 - Check if tumor in center region
 - If yes: drop it
 - If no: extract this patch



Multi Scale Patch Extraction




Data Augmentation

- Augmentations from the paper:
 - Rotate the patch by 0° , 90° , 180° , 270°
 - Apply a left-right flip and repeat rotations
 - Perturb color
 - Small offset of some pixels
- Additional augmentations we did:
 - Apply an up-down flip, since pathology slides do not have canonical orientations


Model Architectures

-


Model Architecture



Single-scale



Multi-scale




Heatmap Construction

Heatmap Construction Methodology

For each testing slide:

- Sliding a window of size 299×299 through the entire image to extract patches
- Using stride = 128 to match the center region's size, so that the prediction do not overlap
- Predict only If the patch contains tissue in its center 128×128 region



Prediction of Patch

For each patch, we calculate prediction result in two ways:


- Method 1:

Do a single Prediction on the patch

- Method 2:

Apply the rotations and left-right flip to obtain predictions for each of the 8 orientations, and average the 8 predictions.

8 Orientations:



Average the 8 predictions

Model Comparison

Single-Scale Models Comparison

Scale Level	Model	Parameters	Validation Accuracy
Level 1	InceptionV3	Transferred	0.9594
		Fine Tuned	0.9891
Level 2		Transferred	0.9297
		Fine Tuned	0.9812
Level 3		Transferred	0.9484
		Fine Tuned	0.9703
Level 1	ResNet50	Transferred	0.80
		Fine Tuned	0.82

Best:
Fine Tuned InceptionV3

Worst:
ResNet50

Thus, we used InceptionV3
model for the following
analysis

Multi-Scale Models Comparison

Scale Level	Model	Parameters	Validation Accuracy
Level 1 & 2	InceptionV3	Transferred	0.9869
		Fine Tuned	0.9906
Level 2 & 3		Transferred	0.9731
		Fine Tuned	0.9859
Level 3 & 4		Transferred	0.9650
		Fine Tuned	0.9859
Level 1 & 3		Transferred	0.9516
		Fine Tuned	0.9641



Best:
Fine Tuned InceptionV3
Using Scale Levels 1 & 2

Worst:
Transferred InceptionV3
Using Levels 1 & 3

We used **finetuned InceptionV3 model** for the following analysis

Model Training Process -- For the best model







- Fine-tuned Interception based model
- Take level 1 and level 2 as input. Level 1 is the reference level that we label the patch.
- We used **early stopping** to prevent overfitting
- **Learning rate** was set to 0.0001



Comparison of Results




This section uses one of the
test data, Slide 091, for
comparison.

Comparison 1: Single-scale v.s. Multi-scale






Comparison 2:


Low Levels v.s. High Levels Scales



Comparison 3: Large v.s. Small Surrounding Context




Metrics




Comparison 4:


Single Prediction v.s. Mean of 8 Predictions Per Patch




True Label




Single-scale Level 1 - Single Prediction
Precision: 0.65, Recall: 0.72
F1 score: 0.69, AUC score: 0.96



Multi-scale Levels 1 & 2 - Single Prediction
Precision: 0.63, Recall: 0.87
F1 score: 0.73, AUC score: 0.97



Single-scale Level 1 - Mean of 8 Predictions
Precision: 0.74, Recall: 0.74
F1 score: 0.74, AUC score: 0.97



Multi-scale Level 1&2 - Mean of 8 Predictions
Precision: 0.60, Recall: 0.58
F1 score: 0.59, AUC score: 0.96

Summary of the Result Analysis

- Multi-scale is better than Single-Scale
- Lower Level is better than Higher Level
- Smaller surrounding area is better than larger surrounding area
- Mean prediction of 8 variations is better than a single prediction of a patch


Best Heatmap Predicting Solution:

- Using **multi-scale model of level 1&2**
- Calculating the **mean prediction of 8 variations** of each patch


Final Prediction

Predicting the heatmaps for
the 3 test slides

Final Predicted Heatmaps on Slide 091




True label




Predicted using Multi-scale
model with level 1&2


Precision score: 0.63
Recall score: 0.87
F1 score: 0.73
AUC score: 0.97



Final Predicted Heatmaps on Slide 110



True label




Predicted using Multi-scale
model with level 1&2

Precision score: 0.96


Recall score: 0.96

F1 score: 0.96


AUC score: 0.99



Final Predicted Heatmaps on Slide 096




True label



Predicted using Multi-scale
model with level 1&2

Precision score: 0.38
Recall score: 0.73
F1 score: 0.50
AUC score: 0.94



**Thanks for
Watching!**



Logic

- Project flow introduction: introduce the overall steps as below
- Data Processing
 - Training & Validation & Testing
 - Some slides has only a little tumor, we dropped them
 - 3 for testing 8 for training and validation
 - Patch extraction for training & validation data
- Model Architectures
 - Transferred
 - InceptionV3 unfix (fine tuned)
 - InceptionV3 unfix (fixed)
 - ResNet unfix (fine tuned)
 - Scales
 - Single Scale: level 1,2,3
 - Multi-scale: level 1&2, 2&3, 3&4, 1&3
- Heatmap Construction
 - Sliding window to extract patches that contain tissue in center part
 - Each patch predict 1 time v.s. predict 8 times?
- Model comparison (表格)
 - Val acc (excel 前两个sheet)
 - Single-scale
 - Multi-scale
 - 结论: 用inceptionV3, 从这之后不show ResNet
- Model comparison w.r.t. heatmaps of 091 (next page)
- Final heatmaps for 3 testing data (091, 110, 096) using the best model

Model comparison w.r.t. heatmaps of 091

- Single or Multi (3 列点点, 横轴是ref level)
 - Inception level 1 v.s. Inception level 1&2
 - Inception level 2 v.s. Inception level 2&3
 - Inception level 3 v.s. Inception level 3&4 不要了
 - 结论: multi的比single好
- Low or high ref level
 - Inception level 1&2 v.s. 2&3 v.s. 3&4
 - 结论: low ref level更好
- How large surrounding area to consider
 - Inception level 1&2 v.s. 1&3
 - 结论: adjacent的更好
- Predict 1 time or 8 times?
 - Inception level 1: 1 time v.s. 8 times
 - Inception level 1&2: 1 time v.s. 8 times
 - 结论: 八次的更好
- Final model:
 - Inception level 1&2 unfix 8 times